

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for DOSYWEL

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and email this form to us on **sales@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



Space Requirements

The weight of the equipment is 3.8kg.

Physical Dimensions

Width	Depth	Height
<i>mm</i>	<i>mm</i>	<i>mm</i>
120	240	220

Environmental Conditions

Altitude: up to 2,000m
 Temperature: from 0 to 45°C
 Relative humidity (RH): from 10 to 80%
 Machine supply voltage: 24 V DC, 3.75 A

Mains Requirements

Electricity Supply

100-240 CAV, 50/60Hz, 2 A

Power Supply: 90W

Socket

1 x Three Pin, 13 amp

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature

Title

Print Name

Establishment